

Silagum

The right solution for all impression techniques



THE RIGHT SOLUTION FOR ALL IMPRESSION TECHNIQUES.

Silagum has been known for years as a proven and reliable impression material based on VPS.

Each single product of the Silagum line is characterized by its exceptional material properties and ease of use.

Silagum's hydrophilicity and thus improved wettability yields a pinsharp reproduction of subgingival areas.

From Handmix over automatically mixing to fully automatic mixing with the MixStar system, products are available in a wide variety of application forms. Silagum offers the right solution for any technique and use.

Good reasons

- Impressions can be stored for a long time
- Excellent recovery after deformation
- Good hydrophilicity
- High-contrast shades
- Tailored to all dental requirements



Excellent recovery after deformation

Due to its high elasticity, Silagum exhibits a very good recovery after deformation, i.e. an extremely low compression set.


It is also this high elasticity that makes removal from the mouth so easy. Any deformation that occurs upon removal safely recovers, thus ensuring precise dimensional reproduction.

Very good hydrophilicity

Through its high surface wettability, Silagum provides for an optimal adaptation to dental hard substances and the mucous membrane in a wet oral environment.

Combined with the flowable viscosity, subgingival areas are precisely reproduced resulting in pin-sharp impressions, even of finest surface details*.

Compression Set in %



- S Silagum-Putty Standard (DMG) P Panasil Putty Soft (Kettenbach)
- H Hydrorise Putty Fast Set (Zhermack)

Source: internal measurement, Data on File

Silagum	Light	Light			Medium			Putty			
	Î		>	Î		Î	>				
		Fast			Fast			MixStar Putty Soft	Standard	Fast	Soft
Mixing time [min]			0:30				0:30		0:30	0:30	0:30
Working time [min]	2:15	1:30	2:15	2:15	1:30	1:30	2:30	2:00	1:45	1:15	1:45
Setting time in the mouth [min]	3:30	2:15	3:30	3:30	2:15	2:30	3:30	3:30	3:30	2:30	3:30
Dimensional change 24 h [%]	≤ 0.4	≤ 0.4	≤ 0.4	≤ 0.5	≤ 0.5	≤ 0.5	≤ 0.35	≤ 0.5	≤ 0.5	≤ 0.5	≤ 0.5
Compression set [%]	≤ 0.4	≤ 0.4	≤ 0.4	≤ 0.5	≤ 0.5	≤ 0.5	≤ 0.4	≤ 0.5	≤ 0.7	≤ 0.7	≤ 0.7
Shore-A-hardness	≈ 50	≈ 50	≈ 50	≈ 50	≈ 50	≈ 60	≈ 54	≈ 62	≈ 70	≈ 70	≈ 62
Strain in compression [%]	≈ 3.5	≈ 3.5	≈ 3.5	≈ 3-5	≈ 3-5	≤ 3.5	≤ 4	1.7-2.6	≈ 1.7	≈ 1.7	≈ 2.5
Color											

Silagum-Putty

Pre-impression material based on VPS in the MixStar cartridge and in the traditional putty jars.

The final hardness of Silagum-MixStar Putty Soft in the MixStar cartridge is adjusted for both an easy removal from the mouth, even with the use of Rimlock trays, and reliable application of correction techniques.

The kneadable pre-impression material Silagum-Putty in the traditional jar is available in three different types.

Good reasons

- Exceptional dimensional
- Snap-set behavior
- Easy and mess-free handling
- Excellent recovery after deformation



Pre-impression with Silagum-Putty.



Pre-impression with relief channels (Silagum-Putty).


Silagum-Light

Proven correction material based on VPS.


The use of Silagum-Light as the second phase of the correction technique or the double mixing technique ensures a precise impression of the marginal area of the preparation. Silagum-Light exhibits a convenient viscosity that simplifies the application, but at the same time ensures safe handling by the user.

Good reasons

- Good hydrophilicity
- Flowable viscosity
- · Available in two versions: regular and fast



*Precise correction impression with Silagum-Putty and Silagum-Light.



*Double mixing impression with Silagum-Putty Soft and Silagum-Light.

Silagum-Mono

Monophase impression material based on VPS.

Silagum-Mono is the material of choice for one-phase and one-step impressions. It is extremely suitable for fixation impressions and also for impressions of crown and bridge preparations.

Good reasons

- Reproduces the finest details owing to excellent flow properties
- Easy and precise application with Intraoral-Tips
- Remains firm in the tray and on the preparation
- Tube material has slight mint flavor



Precise impression with Silagum-Mono.

Silagum-Medium

Medium viscosity impression material based on VPS.

Owing to its consistency of medium viscosity, Silagum-Medium can be used in many ways. It is equally suitable as a correction impression material for two-phase impressions and as the sole material for relining impressions.

Good reasons

- Excellent flow properties
- Available in two versions: regular and fast set



Syringing the preparation with Silagum-Medium for the correction technique.

Silagum-Putty

Package

1 Cartridge @ 380 ml REF 909534 10 MixStar-Tips

2 Jars @ 262 ml base- and catalysatorpaste Putty Standard REF 909018 REF 909037 **Putty Fast Putty Soft** REF 909019

Economy Pack

5 Cartridges @ 380 ml

50 MixStar-Tips REF 909535

8 Jars @ 262 ml

Putty Standard REF 909405 REF 909038 **Putty Fast** REF 909036 Putty Soft

Bulk Pack

16 Jars @ 262 ml

Putty Standard REF 909420 REF 909485 **Putty Soft**

Silagum-Light

Package

2 Cartridges @ 50 ml

12 Automix-Tips, 12 Intraoral-Tips REF 909713 Light Light Fast REF 909714

Bulk Pack

8 Cartridges @ 50 ml

50 Automix-Tips, 50 Intraoral-Tips REF 909715 Light Light Fast REF 909720

8 Tubes @ 90 ml, 2 mixing pads

2 Spatulas REF 909422

Silagum-Mono

Package

2 Cartridges @ 50 ml

4 Automix-Tips

4 Intraoral-Tips REF 909718

Economy Pack

8 Cartridges @ 50 ml 16 Automix-Tips

16 Intraoral-Tips REF 909719

8 Tubes @ 135 ml

2 Mixing pads

2 Spatulas REF 909427

Silagum-Medium

Package

2 Cartridges @ 50 ml

12 Automix-Tips

12 Intraoral-Tips

REF 909716 Medium Medium Fast REF 909717

Accessories

1 Automix-Dispenser Type 50 1:1

REF 909507

50 Automix-Tips Yellow 1:1

REF 909200 REF 909201 50 Intraoral-Tips

1 Bottle @ 10 ml Tray Adhesive

REF 909394

50 MixStar-Tips REF 909567









DMG

Chemisch-Pharmazeutische Fabrik GmbH

Elbgaustraße 248 22547 Hamburg Germany

Fon: +49. (0) 40. 84 006-0 Fax: +49. (0) 40. 84 006-222 info@dmg-dental.com www.dmg-dental.com

