# Luxatemp

## Fluorescence

## Inhalt | Contents

Gebrauchsinformation - DEUTSCH	4
Instructions for use - ENGLISH	14
Mode d'emploi - FRANÇAIS	23
Istruzioni d'uso - ITALIANO	33
Instrucciones de empleo - ESPAÑOL	43
Instruções de uso - PORTUGUÊS	53
Gebruiksaanwijzing - NEDERLANDS	63
Brugsanvisning - DANSK	73
Användarinstruktioner - SVENSKA	82
Instrukcja użycia - POLSKI	91
Информация по применению - РУССКИЙ	101
使用说明书 - 中文	111
<b>사용 정보 - 한국어</b>	<b>119</b>
使用説明書 - 日本語	128

## Produktbeschreibung

---

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## Product Description

---

Luxatemp-Fluorescence is a self-curing composite for the fabrication of temporary crowns and bridges, inlays, onlays and veneers. The 2-component material is based on multifunctional methacrylates and it is automatically mixed. The additional fluorescence makes it possible to create especially natural temporary prostheses with a high esthetic quality.

## Indications

---

- Fabrication of temporary crowns, bridges, inlays, onlays, partial crowns and veneers.
- Creation of long-term temporary prostheses.

## Contraindications

---

Do not use the material in cases of allergies to any of the ingredients or, as can occur in rare cases, in the event of contact allergies.

## Basic safety instructions

---

- For dental use only!
- Keep out of reach of children!
- Avoid contact with skin and eyes! In the event of accidental contact, rinse thoroughly with plenty of water immediately and consult a physician if necessary.
- Wear mouth protection and goggles when constructing the temporary prosthesis, otherwise sanding dust may get into the eyes and respiratory tract!

## Application system

---

- Automix cartridge: see “Using the Automix cartridge” on page 18.
- Smartmix cartridge: see “Using the Smartmix cartridge” on page 20.

**Note:** When using a cartridge or injector for the first time, squeeze out a pea-sized amount of material and discard it.

## Timeline

---

0:00–00:45 min	Insert into the mouth
2:00–3:00 min	Remove from the mouth
6:00–7:00 min	End of curing Work on the temporary prosthesis

**Note:** The times apply at a room temperature of 23 °C and a normal relative air humidity of

50 %. Higher temperatures shorten, lower temperatures increase these times.

## Recommended use

---

1. Before preparing the crown or bridge or prior to extraction of a tooth, take an impression using alginate, silicone or a thermoplastic impression material.

Alternatively, the temporary prosthesis can be made using a laboratory-prepared vacuum-formed stent or, in the case of single crowns, with a strip crown. Adhere precisely to the relevant manufacturer’s instructions!

**Note:** in a silicone casting, undercuts will have to be evened out and if necessary drainage grooves may have to be cut in. To create a connection with the abutment teeth in the case of a gap in the molar region, cut a groove into the casting.

2. Dry the prepared teeth and even out any undercuts in the preparation. Then lightly lubricate the cores and surrounding tissue plus any composite core build-ups using petroleum jelly or a similar separating medium.
3. Select the desired color of Luxatemp-Fluorescence and insert the appropriate cartridge in the Automix-Dispenser or use the corresponding Smartmix injector (see “Using the Automix cartridge” on page 18/ see “Using the Smartmix injector” on page 20).

**Note:** to avoid bubbles, it is important that the tip of the mixing cannula always remains immersed in the material during application.

4. Dispense Luxatemp-Fluorescence into the impression or vacuum-formed stent. Then dispense the material onto the occlusal surfaces of the impression or vacuum-formed

stent and fill up to the gingival areas, only overbuilding slightly.

5. Place the impression or vacuum-formed stent filled with Luxatemp-Fluorescence on the prepared teeth within at most 45s after the start of mixing and hold in place with moderate pressure.

**Note:** the setting reaction has to be monitored intra-orally (e.g. using a probe), because the curing is significantly affected by the temperature in the mouth and the temporary prosthesis can only be removed during its elastic phase.

6. Remove the temporary prosthesis from the mouth during the elastic phase, e. g. 2:00 to 3:00 minutes after the start of mixing.

**Note:** to improve the accuracy of large provisional prostheses further, they can be put back into the impression after removal from the mouth to allow hardening to complete there.

7. Before working on it, remove the oxidation prevention layer from the surface of the temporary prosthesis (e.g. with alcohol).

**Caution!** When processing the temporary, sanding dust can get into the eyes and respiratory tract.

Wear mouth protection and goggles.

Do not inhale the sanded dust and use an extractor.

8. From about 6:00 to 7:00 minutes after mixing, the temporary prosthesis can be processed using cross-toothed grinders, flexible disks, etc.
9. Check occlusion and grind if necessary using suitable instruments.
10. If necessary, the surface of the temporary prosthesis may be glazed with a light-curing high-gloss varnish (e.g. Luxatemp-Glaze & Bond / LuxaGlaze) or polished manually with a rubber polisher or goat's hair wheel.

11. Clean the core, dry it and thoroughly remove the separating medium.

12. Fix the temporary prosthesis.

**Note:** any temporary cement can be using for fixing the temporary. DMG recommends using automatically mixing attachment cements, e.g. TempoCemNE (non-eugenol zinc oxide cement).

Cements containing eugenol can inhibit the polymerization of resin-based attachment cements during the permanent cementation.

13. Check for occlusion and adjust accordingly.

## **Personalization and repair of temporary prosthesis**

---

**Note:** when personalizing and repairing a freshly-made temporary prosthesis it does not need to be roughened and a bonding agent is not required (continue with step 2).

**Caution!** When processing the temporary prosthesis, sanding dust can get into the eyes and respiratory tract.


Wear mouth protection and goggles.  
Do not inhale the sanded dust and use an extractor.

1. Repairing an emplaced temporary prosthesis: Roughen the temporary prosthesis slightly using a grinder or a sandblaster at the appropriate points. Then apply a bonding agent (e.g. Luxatemp-Glaze & Bond / Luxa-Glaze), in accordance with the manufacturer's instructions.
2. Apply Luxatemp-Fluorescence to all surfaces.
3. Fix parts for 3:00 minutes.
4. After 6:00 or 7:00 minutes, the appropriate points can be finished as usual.

**Note:** to speed up the polymerization, the processed temporary prosthesis can be placed in warm water at 50 °C (120°F) for a few minutes.

Alternatively, the temporary prosthesis can also be personalized or repaired using the light-curing, easy-flowing composite LuxaFlow or Luxa-Flow. Follow the associated instructions for use!

## Using the Automix cartridge

### Inserting the cartridge

1. Push up the lever [A] at the back of the Automix-Dispenser and pull back the slide [B] completely.

2. Lift up the plastic lever [C], insert the cartridge and lock it using the plastic lever [C].

**Note:** Make sure that the notches on the cartridge line up with those on the Automix-Dispenser.

### Attaching mixing tip

1. Remove the cartridge cap [D] or used mixing tip by turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

**Note:** Make sure that the notches on the mixing tip line up with those on the cartridge.

3. Turn the mixing tip clockwise 90° to lock it.

### Applying material

When using a cartridge for the first time, squeeze out a pea-sized amount of material and discard it.


The material is mixed in the mixing tip by pulling the trigger [F] of the Automix-Dispenser and it is then applied directly.

**Note:** After the application, leave the used mixing tip on the cartridge to cap unused material!


### Removing the cartridge

1. Push up the lever [A] at the back of the Automix-Dispenser and pull back the slide [B] completely.
2. Lift up the plastic lever [C] and remove the cartridge.

## Using the Smartmix syringe


1. Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it.

2. Attach a new mixing tip [B].

**Note:** Make sure that the notches on the mixing tip line up with those on the Smartmix injector [C].

3. Turn the mixing tip clockwise 90° to lock it.
4. The material is mixed in the cannula when it is expressed and it can be applied directly.

### **Note:**

When using an injector for the first time, squeeze out a pea-sized amount of material and discard it!

Leave the used mixing cannula on the Smartmix injector after use. It serves as a cap.

## Disinfection of the Automix-Dispenser

---

The Automix-Dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfecting solutions listed by the RKI (Robert Koch Institute). The Automix-Dispenser is also autoclavable.

## Side effects

---

There have been no systemic side effects to date.

## Composition

---

Glass filler materials in a matrix of multifunctional methacrylates, catalysts, stabilizers, additives. Free from methyl methacrylate and peroxides. Filler content: 43 wt.% = 24 vol.-%. (0.02 to 1.5  $\mu\text{m}$ )

## Storage and shelf life

---

Store in a dry place at room temperature (15 to 25 °C/59 to 77 °F) !

Do not store in the refrigerator!

Protect against exposure to direct sunlight!

Leave the used mixing cannula on the Smartmix injector after use. It serves as a cap.

Do not use after the expiration date!

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**DMG** Chemisch-Pharmazeutische Fabrik GmbH