# Fill-Up!™

<b>DE</b>	Gebrauchsinformation	3	<b>PL</b>	Instrukcja użycia	68
<b>EN</b>	Instructions for Use	8	<b>SL</b>	Navodila za uporabo	73
<b>FR</b>	Mode d'emploi	13	<b>LT</b>	Naudojimo instrukcija	78
<b>IT</b>	Istruzioni per l'uso	18	<b>ET</b>	Kasutusjuhend	83
<b>ES</b>	Instrucciones de uso	23	<b>CS</b>	Návod k použití	88
<b>PT</b>	Instruções de Utilização	28	<b>HU</b>	Használati utasítás	93
<b>NL</b>	Gebruiksaanwijzing	33	<b>SK</b>	Návod na použitie	98
<b>DA</b>	Brugsanvisning	38	<b>BG</b>	Инструкции за употреба	103
<b>SV</b>	Bruksanvisning	43	<b>LV</b>	Lietošanas instrukcija	108
<b>NO</b>	Bruksanvisning	48	<b>TR</b>	Kullanma Talimatları	113
<b>FI</b>	Käyttöohje	53	<b>RU</b>	Инструкция по применению	118
<b>EL</b>	Οδηγίες χρήσης	58	<b>ZH</b>	使用说明（简体中文）	124
<b>RO</b>	Instrucțiuni de utilizare	63			

**Definition**

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07-2014

## Definition

Fill-Up! is a radiopaque, flowable composite material with antibacterial zinc oxide in an automix syringe and is used for bulk filling therapy in the posterior dental region.

Owing to its dual curing properties, fillings can be administered and cured in a single step in arbitrary thickness.

Etchant Gel S is indicated for etching dentine and enamel for preparatory work for adhesive techniques.

ParaBond Adhesive is a dual-curing, self-conditioning adhesive system for enamel and dentine. It consists of a Non-Rinse Conditioner and a chemically curing adhesive (Adhesives A and B).

## Shades

Fill-Up! is available in the Universal shade. The final colour shade is obtained after 24h and is comparable to a VITA™ A2 / A3.

## Composition

Fill-Up! contains:

Dental glass  
Methacrylates  
Amorphous silica  
Zinc oxide

Etchant Gel S contains:

Phosphoric acid

ParaBond Non-Rinse Conditioner (NRC) contains:

Water  
Acrylamidosulfonic acid  
Methacrylate

ParaBond Adhesive A contains:

Methacrylates  
Polyalkenoates  
Initiators

ParaBond Adhesive B contains:

Ethanol  
Water  
Initiators

## Technical data According to ISO 4049:2009

Average filler particle diameter:	2 µm
Filler particle distribution:	0.1-5.0 µm
Volume filling ratio approx.:	49 %
Weight filling ratio approx.:	65 %

## Indications

- Direct restoration of Class I and II cavities
- Cavity lining - first layer for Class I and II cavities
- Core build-ups

## Contra-indications

Proven hypersensitivities to ingredients contained in Fill-Up! or in the mentioned complementary products. Inadequate oral hygiene. If it is impossible to keep the working area dry.

## Side effects

Fill-Up! ingredients can cause sensitivity in predisposed persons.

## Interactions with other agents

Phenolic and other substances (e.g. zinc oxide eugenol) which inhibit polymerisation must not be allowed to come into contact with Fill-Up!.

If brought into contact with the etched and dried enamel, saliva or other liquids will destroy the adhesive system's bonding effect.

## Clinically measured times

Polymerisation begins upon first contact between base and catalyst.

	Intraoral 37 °C / 99 °F
Working time	1 : 00 min
Setting time (incl. working time)	3 : 00 min

The processing time depends on the temperature. This decreases at higher temperatures and is prolonged at lower temperatures. Fill-Up! is light-sensitive and should not be exposed to intense light for longer

than 30 s prior to curing, particularly to surgical light.

## Application

Cleansing the tooth

Clean the tooth to be restored and the mesial and distal adjacent teeth with a brush and fluoride-free cleaning paste.

Tooth isolation

A dry working area is the prerequisite for outstanding results. The use of a rubber dam is recommended (e.g. HYGENIC® or ROEKO Dental Dam).

Preparing the cavity

When preparing the cavity one should employ techniques gentle to the hard tooth substance wherever possible (principle of „adhesive restoration technique“).

Use 80 µm preparatory diamonds for working on enamel and dentine and finish with 25 µm diamonds. Beveling of the enamel margin is recommended to increase the bonding area between the tooth and the filling material, to prevent fractures of the enamel margin and thus to optimise the marginal fit.

Application of Etchant Gel S

1. Directly apply Etchant Gel S to dentine and enamel with the applicator needle.
2. Etch for 15 s and afterwards rinse thoroughly with water until all traces of Etchant Gel S have been removed.
3. Dry with pressurised oil-free air, taking care not to overdry the dentine. Etched enamel has a chalky white appearance.

*Note: Gloves, safety goggles (for both doctor and patient) and mouth guards are important for safety reasons! Avoid bringing Etchant Gel S into contact with gingiva or adjacent teeth. Use fluoride on enamel which has been etched and is not covered by composite.*

*As an alternative to etching with Etchant Gel S, the cavity can be treated with ParaBond Non-Rinse Conditioner for 30 s and dried with oil-free air for 2 s.*

#### Application of ParaBond Adhesive

1. Mix a drop of Adhesive A with a drop of Adhesive B.
2. Massage the mixed adhesive into the cavity and on the dentine and enamel with a brush for 30 s.
3. Remove excess adhesive with a gentle stream of air for 2 s.

**Important:** The time between applying the adhesive and filling the cavity must not be longer than 5 min. Otherwise repeat the „application of ParaBond adhesive“.

#### Matrix / Interdental wedges

Apply a thin, cambered matrix band. Fixate the matrix approximately with interdental wedges.

## **A. Direct restoration of Class I and II cavities**

### **Filling without a covering layer**

#### Application of Fill-Up!

1. Remove the safety cap or mixing tip. Press a small amount of the material onto a paper towel until the base and catalyst both flow in equal amounts from the opening. This ensures a homogenous mixture.

2. Wipe the opening of the syringe with a paper towel. Attach the mixing tip and tighten with ¼ turn clockwise (90°). Squeeze material and discard until even, homogenous paste flows.
3. Place mixing tip in the deepest section of the cavity and stay immersed in the material during application to avoid air pockets. Fill the entire cavity in a single step.

**Important:** to avoid air pockets, do not shape during the processing period. The surface of the cavity must be completely cured before using rotating instruments.

4. After use, rub mixing tip with disinfectant and do not remove.

#### Polymerisation

Curing commences with the first contact between base and catalyst. After 3 min, Fill-Up! is chemically cured.

To accelerate the curing process and obtaining a surface suitable for grinding, curing by light can be performed as follows:

Light output of the polymerisation lamp	Exposure time
800 mW/cm <sup>2</sup>	10 s
1200 mW/cm <sup>2</sup>	7 s
1600 mW/cm <sup>2</sup>	5 s

*Note: the recommended exposure times should not be exceeded to minimise shrinking stress. It is also recommended to comply with the curing time in the case of very deep cavities to avoid bonding of the material to the matrix.*

#### Finishing

Finishing includes three steps: rough work-up, fine work-up and polishing. Sets of finishing diamonds (40 µm and 15 µm) can be used universally and are the least destructive. Application should be without any pressure, with a constant wiping motion and copious spraying with water. The ideal speed lies between 5000 and 15 000 rpm. Flexible, diamonded files (80 µm, 40 µm and 15 µm) and strips are suitable for processing approximal surfaces. To obtain perfect polishing results, use, for example, the DIATECH composite polishing instruments or the SwissFlex polishing system. It is advised to fluoride tooth substance which is not covered by composite.

### **Filling with a covering layer**

#### Application of Fill-Up!

Proceed as described under „Filling without a covering layer - application of Fill-Up!“.

**Important:** do not apply Fill-Up! up to the occlusal preparation margin in order to leave space for the covering layer.

To reduce the inhibition layer of Fill-Up!, subsequent light curing is recommended (see „Filling without covering layer - polymerisation“).

#### Application of a covering layer

**Important:** the inhibition layer must not be contaminated. Should this occur, repeat the application of ParaBond Adhesive A + B.

Place the covering layer with a methacrylate-based universal/posterior tooth composite (preferably composites by Coltène/Whaledent AG) and light-cure according to the manufacturer's instructions.

#### Finishing

Proceed as described under „Filling without a covering layer - finishing“ or according to the manufacturer's instructions.

## **B. Cavity lining**

Fill-Up! can be used as a cavity liner similar to a flowable composite.

## **C. Core build-up**

#### Application of Fill-Up!

1. Remove the safety cap or mixing tip. Press a small amount of the material onto a paper towel until the base and catalyst both flow in equal amounts from the opening. This ensures a homogenous mixture.
2. Wipe the opening of the syringe with a paper towel. Attach the mixing tip and tighten with ¼ turn clockwise (90°). Squeeze material and discard until even, homogenous paste flows.
3. If required, place a matrix strip around the prepared tooth.
4. Apply Fill-Up! directly onto the preparation. After use, wipe off the mixing tip with disinfectant and do not remove.

*Note: it is important for adequate healthy tooth substance to be present, to achieve a suitable ferrule effect (1.5 mm apical around the stump).*

5. To accelerate the curing process and/ or to reduce the inhibition layer, light can be used to cure for 20 s.

## **Emergency measures**

If the product comes into contact with the oral mucosa, simply rinse

with water. If the product gets in the eyes, wash them out thoroughly with water (10 min) and then consult an ophthalmologist.

### Notes

Only to be provided to dentists and dental technicians or on their behalf. Keep out of the reach of children. Do not use the product beyond its expiry date. All bottles must be sealed immediately after each application.

### Shelf life and labelling

The expiry date and the LOT identification can be found on the containers and the outer packaging.

### Storage

Store Fill-Up! and ParaBond Adhesive in the refrigerator (4 – 8°C / 39 – 46°F). After the first application and in case of timely consumption, storage is possible at room temperature (approx. 23°C / 73°F). Store Etchant Gel S at room temperature (approx. 23°C / 73°F). Do not expose materials to direct sunlight or other heat sources.

### Caution

Federal law restricts this device to sale by or on the order of a dentist.

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**Coltene/Whaledent AG**