# Venus® Pearl

- GE** Gebrauchsanweisung
- IB** Instructions for use
- RD** Mode d'emploi
- ES** Instrucciones de uso
- FI** Istruzioni per l'uso
- FR** Instruções de uso

6031885/00

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## Venus® Pearl

**Universal light curing nano hybrid composite, radio-opaque**
**Venus Pearl** is a light curing, radio-opaque nano composite. It is indicated for adhesively bonded, tooth coloured restorations. It is suitable for anterior and posterior applications. Due to the excellent material properties, a beautiful aesthetic result and a durable, high lustre polish is easily reached. Its radiopacity is equivalent to 2.0 - 2.5 mm of aluminium (dentine = 1 mm, enamel = 2 mm).

**Composition**
Nanofüllstoffe: Bio-Flur-Silicatglas, TiO2-Urethansäurelylate, Silica, UDMA, TEGDMA, polymer, Triäthoxyd, fluoreszierende Pigmente, Metallische Oxide Pigmente, organische Pigmente, Antioxidanzien, BISF, Kampheraromen. Venus Pearl contains approximately 59% fill by volume, with 50% organic fill by volume. Its fill particle size is 0.5 - 5 µm.

**Barium Aluminium Bio-Flur Silicat Glas**, TiO2-Urethansäurelylate, Silica, UDMA, TEGDMA, polymer, Triäthoxyd, fluoreszierende Pigmente, Metallische Oxide Pigmente, organische Pigmente, Antioxidanzien, BISF, Kampheraromen. Venus Pearl contains approximately 59% fill by volume, with 50% organic fill by volume. Its fill particle size is 0.5 - 5 µm.

• Direct restoration of Class I – V cavities (including class IV, VI, Black)
• Direct composites veneers
• Shape corrections of teeth (i.e. diastemas, congenital defects in teeth, etc.)
• Splinting of loose teeth
• Closing traumatic or periodontal bed damaged margins
• Indirect restorations (inlays, veneers)
• Restoration of primary teeth
• Core build-up in repairs of porcelain, composite (in combination with an adequate repair-system)

**Kontraindikationen:**
Application on pulp exposed dentin (less than 1 mm).
Use of this product is contraindicated in case of known or suspected allergies against components of this product. Do not use this product in case of known or suspected allergies against methacrylate compounds.

**Side effects**
This is not one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

**Warnings and safety instructions:**
• May be irritating to eyes and skin. Avoid contact with skin. Use caution when sensitization by skin contact. Wear suitable protective goggles. In case of contact with skin, immediately wash with water and mild soap. Avoid contact with eyes. In case of contact with eyes, rinse with water for several minutes. If symptoms persist, consult an eye specialist.

**Application hints:**
• Visually inspect for damage before use. Damaged products must not be used. Composites should only be used at room temperature.

• Do not use deep cavity cavities. Use only for minimal penetration of pulp which may be clearly clinically evident i.e. where the pulp fills a small cavity (drill only through the pulp) or where the pulp is exposed to the cavity (drill only through the pulp) or where the pulp is exposed to the cavity (drill only through the pulp).

• Do not use beyond contact margin (around the cavity preparation) because there are risks of overexposure. Do not use after composite finishing.
• Use protective eyewear for sprayer protection or protection of patients.
• The PL1 is not intended for single use only or reprocessing or re-use as is not possible without damaging and contaminating the restorative material by those phenolic disinfection solutions. Do not use for protection of remaining restorative material.
• The PL1 is not intended for single use only or reprocessing or re-use as is not possible without damaging and contaminating the restorative material by those phenolic disinfection solutions. Do not use for protection of remaining restorative material.

• When using Venus Pearl PL1, use a gentle, even motion when exerting pressure on the applicator gun. Use of excessive force or sudden movement could result in unanticipated extrusion of material or fracture of the PL1 from the applicator gun, creating a potential hazard. We recommend using the Kutzer PL1 Applicator STEEL.
When using other PL1 applicator guns see full instructions and precautions of these devices.

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