Instruction Manual GenTek™ Hexalobular Screw*

*Incl. GenTek™ Hexalobular Titanium and Gold-Tite® Screws

Version: 3 / 01. 2022
Compiled at: 28. 01. 2022
Effective date: 01. 03. 2022
Replaces version: 2 / 05. 2021

CONTENT

EN (English) Instruction Manual GenTek™ Hexalobular Screw	2
BG (Български) Ръководство за употреба Винт GenTek™ с глава с гнездо „шестолъчна звезда“	6
CS (Česky) Návod k použití GenTek™ šestihranný šroub	10
DA (DANISH) Brugervejledning GenTek™ hexalobulær skrue	14
DE (Deutsch) Gebrauchsanweisung GenTek™ Sechsrundschrauben	18
EL (Ελληνικά) Εγχειρίδιο οδηγιών GenTek™ Εξάλοβο κατσαβίδι	22
ES (Español) Manual de instrucciones GenTek™ Tornillo Hexalobular	26
ET (Eesti) Kasutusjuhend GenTek™ Hexalobular Screw	30
FI (Suomi) Käyttöohje Heksalobulaarinen GenTek™ -ruuvi	34
FR (Français) Mode d'emploi Vis hexalobée GenTek™	38
HR (Hrvatski) Priručnik GenTek™ heksalobularni vijak	42
HU (Magyar) Kezelési útmutató GenTek™ Torx Csavar	46
IT (Italiano) Istruzioni d'uso Vite esalobata GenTek™	50
LT (Lietuvių k) Naudojimo instrukcija GenTek™ šešiabriaunis varžtas	54
LV (Latviski) Lietošanas rokasgrāmata GenTek™ zvaigžņveida skrūve	58
NL (Nederlands) Gebruikershandleiding Hexalobulaire GenTek™ -schroef	62
PL (Polski) Instrukcja stosowania Śruba torx GenTek™	66
PT (Português) Manual de Instruções GenTek™ Parafuso Hexalobular	70
RO (Română) Manual de instrucțiuni Șurub hexalobat GenTek™	74
SK (Slovak/Sloavkian) Návod na použitie Skrutka s hexalobulárnou drážkou GenTek™	78
SL (Slovenščina) Navodila za uporabo Heksalobularni vijak GenTek™	82
SV (Svenska) Bruksanvisning GenTek™ hexalobulär skruv	86

Instruction Manual GenTek™ Hexalobular Screw* *Incl. GenTek™ Hexalobular Titanium and Gold-Tite® Screws

Version: 3/01.2022
Compiled at: 28.01.2022
Effective date: 01.03.2022
Replaces version: 2/05.2021

Indications

Intended Use / Intended Purpose

The GenTek™ Hexalobular Screws are intended to retain GenTek™ Angulated Screw Channel TiBases and custom abutments on screw and cement-retained single-unit (engaging), multi-unit (non-engaging), and full-arch prosthetic restorations to the implant (Final Prosthetic Screw) or analog (Laboratory Try-in Screw).

The GenTek™ Angulated Screw Channel TiBase systems are intended to be used only with the corresponding GenTek™ Hexalobular Screws (final prosthetic screw and Try-in Screw).


Device description

GenTek™ Hexalobular Titanium Screws for TSV®/Trabecular Metal™ and Eztetic® connections are packaged individually and also with the GenTek™ Angulated Screw Channel TiBases. GenTek™ Hexalobular Titanium Screws for TSV® and Eztetic® serve both as the final prosthetic screw and the laboratory Try-in Screw.

GenTek™ Hexalobular Gold-Tite® Screws for Certain® connections are packaged individually and also with the GenTek™ Angulated Screw Channel TiBases. GenTek™ Hexalobular Gold-Tite® Screws for Certain® connections serve as the final prosthetic screw and the GenTek™ Hexalobular Titanium Try-in Screws serve as the laboratory screw.

GenTek™ Hexalobular Gold-Tite® Screws for External Hex connections are packaged individually. GenTek™ Hexalobular Gold-Tite® Screws for External Hex connections serve as the final prosthetic screw for a custom abutment with an external hex connection. The GenTek™ Hexalobular Titanium Try-in Screws serve as the laboratory screw.

The GenTek™ Hexalobular Titanium and Gold-Tite® Screws are for single use only. All the Hexalobular Screws are provided non-sterile and intended to be cleaned and sterilized prior to placement into the patient's mouth.

Hexalobular Drive

GenTek™ Hexalobular Screws allow the placement of the screw channel up to 25°. Hexalobular Screws contain a hexalobular drive that is intended to be used with a corresponding GenTek™ Hexalobular Screwdriver.

Materials

- × GenTek™ Hexalobular Titanium Screws are manufactured from titanium alloy grade 5 ELI (TiAl6V4 ELI)
- × GenTek™ Hexalobular Gold-Tite® Screws are manufactured from stainless steel 316L and contain a 99.9% pure gold surface coating

Procedure

- × Placing/Removing the GenTek™ Screws


GenTek™ Hexalobular Titanium Screws TSV®/Trabecular Metal™ and Eztetic®

Note: GenTek™ Hexalobular Titanium Screws (TSV®/Trabecular Metal™ and Eztetic®) are utilized for both Laboratory (Try-in) and Clinical (Final) settings for both single-unit and multi-unit restorations. The screw used for try-in in the laboratory setting should never be re-used for final patient restoration.


Placing in a Laboratory/Clinical Setting

1. Fully insert the Hexalobular Titanium Screw into the channel of restoration (prior to placing into the analog/implant) using the Hexalobular Screwdriver. To ensure that the screw is properly inserted, the screw threads must be visible at the bottom of the connection.

Note: The screw head external geometry for the GenTek™ Hexalobular Titanium Screws (TSV®/Trabecular Metal™ and Eztetic®) has 3 small protrusions.

2. Once the Hexalobular Screw is fully inserted into the TiBase, keep pressure on the screw with the Hexalobular Screwdriver while placing the TiBase/Screw assembly onto the implant in the desired orientation. Tighten the TiBase/Screw assembly to the implant using the driver. Refer to the below table for placement torque recommendations for Hexalobular Screws.

Instruction Manual

GenTek™ Hexalobular Screw* *Incl. GenTek™ Hexalobular Titanium and Gold-Tite® Screws

Version: 3 / 01. 2022
Compiled at: 28. 01. 2022
Effective date: 01. 03. 2022
Replaces version: 2 / 05. 2021

Removing the GenTek™ Hexalobular Titanium Screws

- Should the restoration/screw (TSV®/Trabecular Metal™ and Eztetic®) system require disassembly from the implant, an abutment removal tool is not required. The Hexalobular Titanium Screw is intended to assist in the removal of the abutment during standard screw removal.

Note: Should the restoration remain assembled to the implant once the screw has been loosened, direct manual manipulation of the restoration with or without a tool such as a hemostat, may be used to facilitate full removal.

Note: During the removal of the TiBase/Screw (TSV®/Trabecular Metal™ and Eztetic®) system the clinician may feel resistance twice. There will be resistance when initially loosening the screw torque and then again when the screw engages the TiBase to separate the system.

GenTek™ Hexalobular Gold-Tite® Screws Certain® and External Hex

Note: There are several GenTek™ Hexalobular Gold-Tite® Screws:

Connection	Restoration Type
Certain®	Engaging / Non-Engaging
External Hex	Engaging

Placing in a Laboratory Setting


Note: GenTek™ Hexalobular Gold-Tite® Screws are for final prosthetic use and are not intended for laboratory use.

Refer to "Placing the GenTek™ Hexalobular Titanium Try-in Screws" section in this IFU.

Placing in a Clinical Setting

Once the Restoration is seated into the implant, place the Hexalobular Gold-Tite® Screw into the channel of the Restoration and tighten the Restoration/Screw assembly to its prosthetic torque using the Hexalobular Screwdriver. Refer to the table for placement torque recommendations for Hexalobular Screws.


GenTek™ Hexalobular Titanium Try-in Screws Certain® and External Hex


Note: There are several GenTek™ Hexalobular Titanium Try-in Screws:

Connection	Restoration Type
Certain®	Engaging / Non-Engaging
External Hex	Engaging

- When placing the laboratory Try-in Screw onto the model, use the corresponding restoration to ensure the proper assembly will be affixed to the analog.

- Fully insert the Hexalobular Titanium Try-in Screw into the channel of the restoration, using the Hexalobular Screwdriver and hand-tighten.

Instruction Manual GenTek™ Hexalobular Screw*

*Incl. GenTek™ Hexalobular Titanium and Gold-Tite® Screws

Version: **3 / 01. 2022**
Compiled at: **28. 01. 2022**
Effective date: **01. 03. 2022**
Replaces version: **2 / 05. 2021**

Note: Only corresponding GenTek™ Hexalobular Screws for each connection (Certain®, External Hex, TSV®/ Trabecular Metal™ and Eztetic®) must be used with the corresponding restorations containing an angulated screw channel. Screws from other GenTek™ systems must not be used with restorations containing an angulated screw channel.

Refer to the below table for placement torque recommendations for GenTek™ Hexalobular (laboratory try-in or final prosthetic) Screws:

Connection*	Angulation Range	Laboratory Try-in Screw	Final Prosthetic Screw Torque	GenTek™ Hexalobular Driver
Certain®	0–25 degrees	Hand-tighten**	20 Ncm	ZFX02HLD21, ZFX02HLD28, ZFX02HLD36
TSV®/ Trabecular Metal™	0–25 degrees	Hand-tighten**	30 Ncm	ZFX02HLD21, ZFX02HLD28, ZFX02HLD36
Eztetic®	0–25 degrees	Hand-tighten**	30 Ncm	ZFX02HLD21, ZFX02HLD28, ZFX02HLD36
External Hex	0–25 degrees	Hand-tighten**	32 Ncm	ZFX02HLD21, ZFX02HLD28, ZFX02HLD36

*Applicable to single-unit and multi-unit for each connection
**Minimum of 5 Ncm; Maximum of 15 Ncm

Cleaning and Sterilization

The GenTek™ Hexalobular Screws are non-sterile. Prior to placement of the GenTek™ Hexalobular Screw into GenTek™ Angulated Screw Channel TiBase / restoration and implant, the GenTek™ Hexalobular Screw must first be cleaned and sterilized.

Step 1: Manual Cleaning

1. Rinse the instruments under running cold tap water for 2 minutes.
2. Hand wash instruments in a bath with demineralized water for 5 minutes. Remove debris with a soft brush. Cavities should be rinsed with a syringe. This step should be repeated until all debris is removed.
3. Clean the instrument in an ultrasonic bath using a mildly alkaline cleaner (1.5%) at a frequency of at least 40 kHz for 20 minutes.
4. Rinse thoroughly with demineralized water for 1 min.

Step 2: Sterilization

One of the following sterilization methods is recommended to be used:

- × **Method 1:**
Full cycle pre-vacuum steam sterilization at a temperature of 134°C (273.2°F) for 3 minutes, 4 pulse minimum.
- × **Method 2: (alternative to method 1)**
Full cycle gravity steam sterilization at 134°C (273.2°F) for 3 minutes minimum.

Step 3: Drying

After sterilization treatment it is recommended to dry the products at a temperature of maximum 134°C (273.2°F) for 10 minutes. The products are ready for use 5 minutes after completion of the drying process.

Note: Follow the loading instructions issued by the sterilizer manufacturer.

Storage and Handling

The GenTek™ Hexalobular Screws should be stored at room temperature.

Contraindications

None

Known Side Effects

Allergies to the alloy or contents of the alloy may very rarely occur.

Warnings

- × The GenTek™ Hexalobular Screws must only be used and handled by dental professionals who have approved training in dentistry.
- × GenTek™ Hexalobular Screws may only be used for Zimmer Biomet Dental implant systems.
- × GenTek™ Hexalobular Screws are small and, therefore, they must be handled with caution to avoid swallowing or inhaling by the patient.
- × The use of an incorrect prosthetic screw may result in damage to the GenTek™ Angulated Screw Channel TiBase, the prosthetic screw and/or the implant.
- × The use of a prosthetic screw torque higher than the recommended value may result in damage to the GenTek™ Angulated Screw Channel TiBase, the GenTek™ Hexalobular Screw and/or the implant. The use of torque values lower than those recommended may result in loosening of the GenTek™ Angulated Screw Channel TiBase and/or the GenTek™ Hexalobular Screw.
- × The GenTek™ Hexalobular Screwdriver must be in good condition. A worn screwdriver can lead to over-rotation of the hexalobe in the screw head which can result in the screw not being fully tightened and/or not being able to be removed.
- × For any lab work and any try-in of the GenTek™ Angulated Screw Channel TiBase, a corresponding lab screw must be used. The included prosthetic screw may only be used to permanently fix the finished restoration.

Instruction Manual GenTek™ Hexalobular Screw*

*Incl. GenTek™ Hexalobular Titanium and Gold-Tite® Screws

Version: 3 / 01. 2022
Compiled at: 28. 01. 2022
Effective date: 01. 03. 2022
Replaces version: 2 / 05. 2021

Explanation of symbols

Non-sterile

Do not use if package is damaged

Manufacturers Information Address of legal manufacturers facility

Don't use twice

Recommended torque to be applied for final placement of the included screw

Recommended to fix the Try-In Screw hand-tight

LOT-Number, Symbol followed by Production Date of the lot

Medical Device

Reference Number, Symbol followed by Manufacturer's Item Number

Quantity / Unit of package content

Unique Device Identifier

Consult instructions for use www.zfx-dental.com

Consult instructions for use for mating hexalobular screws and angulated screw channel TiBase abutments.

Information on Manufacturer

Zfx GmbH
Kopernikusstraße 15
85221 Dachau
Germany
T +49 (0) 8131 / 33 244 - 0
F +49 (0) 8131 / 33 244 - 10
www.zfx-dental.com

Incident Reporting / Complaints

Please send to emeacomplaints@ZimVie.com

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4094DE REV. 01/22 ©2022 Zfx GmbH & Zimmer Dental Inc. | Zimmer und das Zimmer-Logo sind Marken von Zimmer Inc. oder ihrer Tochtergesellschaften und Zfx und das Zfx-Logo sind Marken von Zfx GmbH, Deutschland. Alle anderen Marken sind Eigentum des jeweiligen Besitzers.