NX3

**NEXUS®
THIRD
GENERATION**

**Instructions
For Use**

Kerr™

NX3 Universal Resin Cement System

NX3 is a permanent resin cement system that offers innovative chemistry for unmatched esthetics, delivery system versatility and compatibility with both total-etch and self-etch adhesives. Delivery system choices include an automix dual syringe dual-cure cement (base/catalyst) and a single-syringe light-cure (base only) cement for multiple veneer units where flexible work time is needed. NX3 dual-syringe can be used for all indirect applications including veneers and metal-based restorations, without the need of an activator for the adhesive. A proprietary amine-free initiator system provides excellent color stability for long-term esthetics. One set of try-in gels match both light-cure and dual-cure resin cements. NX3 contains nano filler technology for improved dispensing, handling, cleanup, and polishing properties.

Indications for use:

Cementation of veneers, inlays, onlays, crowns, bridges, and posts.

Adhesive bonding of amalgam restorations.

Core-buildup material.

Restorative materials:

Indicated for cementation of ceramic, porcelain, resin, metal-based materials and CAD/CAM blocks.

Note: When a total-etch technique is desired, the use of OptiBond™ Solo Plus adhesive with NX3 resin cement is recommended. When a self-etch technique is desired, the use of OptiBond XTR*** is recommended.**

** or OptiBond S

*** or OptiBond Versa

I. Fit Evaluation of Restoration, Try-in, and Restoration Surface Preparation

1. Remove temporary restoration and thoroughly clean the tooth (teeth) with fluoride-free prophylactic paste. Rinse thoroughly with water and lightly dry with air.
2. Evaluate fit. Make any necessary adjustments.
3. To evaluate shade (for composite and ceramic restorations), place a thin layer of the selected water-soluble try-in paste on the restoration. **The try-in paste matches both the single-syringe light-cure cement, and the dual-syringe dual-cure cement.** Using light pressure, completely seat the restoration. Remove any excess try-in paste.
4. Evaluate shade. Remove restoration (s). Wash out the water-soluble paste with a strong water spray. Remove any organic debris by cleaning with alcohol or acetone, followed by cleaning in water in an ultrasonic cleaner. The ceramic or porcelain surface may be additionally cleaned with phosphoric acid etchant followed by thorough rinsing and drying.
5. Prepare the internal surface of the restoration according to manufacturer's instructions. Here are some general recommendations:

Ceramic/Porcelain/Composite Restorations: Sandblast the surface with 50m alumina (with a pressure of about 15 psi (0.1 MPa) for composite restoration or about 30 psi (0.2 MPa) for porcelain/ceramic restoration) and etch the restoration with hydrofluoric acid for 1 minute. Apply silane to the internal portions of the restoration (s). Lightly air thin. Place the restorations in a lightproof box to await the cementation procedure. Kerr silane contains resin, so ambient light can cause premature curing. The resin in the silane also precludes the need for the application of additional bonding agent on the restoration prior to cementation. When OptiBond XTR is being used, OptiBond XTR can be used in place of Kerr Silane. Light curing of OptiBond XTR ADHESIVE is optional.

Metal-based, Zirconia-based, Alumina-based restorations: Sandblast the internal surface with 50m alumina with a pressure of about 60 psi (0.4 MPa). Tin plating of high noble or gold metals is not required but may further enhance the adhesion. The metal surface is primed with either OptiBond Solo Plus adhesive or OptiBond XTR ADHESIVE by applying a coat of adhesive on the internal surface of the restoration and air thin the adhesive. Light-curing of adhesive is optional. Place the restoration in a lightproof box to await the cementation procedure.

II. Tooth Preparation

Either OptiBond Solo Plus (**total-etch**) or OptiBond XTR (**self-etch**) adhesives can be used in combination with NX3.

Total-Etch Technique

OptiBond Solo Plus Directions:

1. After thoroughly cleaning the preparations (pumice and prophyl cup), etch enamel and dentin for 15 seconds with Kerr Gel Etchant (37.5% phosphoric acid). Rinse thoroughly and air dry (or blot dry). Do not desiccate.
2. Apply OptiBond Solo Plus to the enamel/dentin surfaces with an applicator tip for 15 seconds using a light brushing motion.
3. Air thin the adhesive for 3 seconds. Avoid pooling of adhesive before light curing (for inlay, onlay, crown, and post cementation). **Remove excess adhesive with a dry applicator brush or absorbent paper point (for post cementation).**
4. Light cure for 10 seconds.*

Self-Etch Technique

OptiBond XTR Directions:

1. Thoroughly clean the preparations (pumice and prophyl cup). Wash thoroughly with water spray and air dry. Do not desiccate.
2. Using the disposable applicator brush, apply PRIMER to the enamel/dentin surface. Scrub the surface with a brushing motion for 20 seconds. Air thin for 5 seconds with medium air pressure.
3. Apply ADHESIVE to the enamel/dentin surface with a light brushing motion for 15 seconds. Air thin with gentle air first and then strong air for at least 5 seconds to avoid pooling of adhesive before light curing.
4. Light-cure for 10 seconds.* (Light-curing of OptiBond XTR ADHESIVE is optional. When light-curing of OptiBond XTR ADHESIVE is omitted, OptiBond XTR ADHESIVE needs to be air thinned with maximum air pressure for 15 seconds prior to NX3 application).

III. Cement Application and Seating of Restoration

Work Time – Self-cure mode, a minimum of 1.5 minutes.

Clean up – Approximately 2-3 minutes after placement, or tack cure for 2 seconds for immediate cleanup.

Note:

1. Work time is based on material at 23°C. Clean-up time is based on material at 37°C. This time may vary based on storage conditions, temperature, humidity, and age of product.
2. For the dual-cure cement, bleed cartridge before initial use.

A. Veneers

Note: Both single-syringe light-cure resin cement and dual-syringe, dual-cure resin cement (ideal for 1-2 units) can be used for cementing the veneers. For thicker veneers, a dual-cure resin cement is recommended.

Dispense cement directly onto the internal surface of the veneer. Place the veneer gently on the tooth, allowing the cement to slowly flow from all margins. Spot cure the veneer into place on the facial surface away from the margins using a small diameter light guide for 10 seconds. After removing the excess cement, light cure all surfaces for a minimum of 20 seconds* per surface.

B. Crowns, bridges, inlays, onlays, and metal-based restorations

Note: On the prep, avoid pooling of adhesive before light curing of adhesive. Excess adhesive can be removed with a dry applicator brush before light-curing of the adhesive.

Apply the dual-cure cement to the restoration or the prep. Seat the restoration gently onto the preparation allowing the cement to flow from all sides. Remove excess cement.† Light cure all surfaces for a minimum of 20 seconds* per surface.

† Excess cement cleanup – Excess cement is best removed in its gel state with a scaler or explorer. Gel state can be achieved by tack curing excess with a light for 1-2 seconds or allowing the cement to self-cure for 2-3 minutes after application.

C. Post Bonding and Core buildup

Note: On the prep, avoid pooling of adhesive before light curing of adhesive. Excess adhesive can be removed with a dry applicator brush or absorbent paper point before light-curing of the adhesive.

1. Prepare post space. Size and fit post.
2. Place adhesive on the post in a thin even layer, lightly air thin if necessary. Light-curing of the adhesive is optional.
3. Apply the dual-cure cement onto the post and/or into the post preparation, seat the post, and vibrate the post slightly to avoid the possibility of trapped air.
4. When you are satisfied that the post is properly seated, remove all excess cement. Light cure all surfaces for a minimum of 20 seconds* per surface.
5. Proceed with core buildup and/or preparation according to manufacturer's instructions for use.

Note: NX3 can be used as a core buildup material.

IV. Finishing and Polishing

1. Remove polymerized excess with a finishing diamond bur.
2. Use finishing strips for proximal area.
3. Polish margins with Kerr Gloss Plus discs, cups, or points. Use Kerr HiLuster Plus Dia Polishers for a final high luster shine.
4. Check occlusion and make adjustments if necessary.

***Recommended Cure Times:** Demi/Demi Plus, 5 seconds; L.E.Demetron II, 5 seconds; L.E.Demetron I, 10 seconds; Optilux 501, 10 seconds. For all other lights, see manufacturer's recommendation.

CAUTION

1. For phosphoric acid gel etchant and adhesive, avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately and thoroughly with water. Get medical attention for eyes. Do not take internally.
2. Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes, and soft tissue. Wash thoroughly with water after contact.

Storage and Shelf Life

Store NX3 at ambient temperature.

NX3 is not recommended for use beyond the stated expiration date marked on the packaging.

"OptiBond" is a trademark owned by Kerr Corporation.

"OptiBond" is registered with the U.S. PTO and other countries.

Limited Warranty - Limitation of Kerr's Liability

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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Kerr Italia, S.r.l.
Kerr Australia Pty. Ltd.
Sybron Dental Specialties Japan, Inc.