## CONTENT

<b>EN (English)</b> Instruction Manual GenTek™ Angulated Screw Channel TiBase*	2
<b>BG (Български)</b> Ръководство за употреба Титанова база GenTek™ TiBase с ангуларен винтов канал	6
<b>CS (Česky)</b> Návod k použití GenTek™ TiBase se zaúhleným šroubovým kanálem	10
<b>DA (DANISH)</b> Brugervejledning GenTek™ TiBase med vinklet skruekanal	14
<b>DE (Deutsch)</b> Gebrauchsanweisung GenTek™ TiBasen mit abgewinkelttem Schraubenkanal	18
<b>EL (Ελληνικά)</b> Εγχειρίδιο οδηγιών Στήριγμα τιτανίου GenTek™ με γωνιώδη αυλό βίδα	23
<b>ES (Español)</b> Manual de instrucciones Base de titanio GenTek™ con canal de tornillo angulado	28
<b>ET (Eesti)</b> Kasutusjuhend GenTek™-i -i anguleeritud kruvikanal TiBase	33
<b>FI (Suomi)</b> Käyttöohje GenTek™ TiBase -kanta, joissa on kallistettu ruuvikanava	37
<b>FR (Français)</b> Mode d'emploi TiBase à canal de vis angulé GenTek™	41
<b>HR (Hrvatski)</b> Priručnik GenTek™ titanska baza s kanalom za vijak s nagibom	45
<b>HU (Magyar)</b> Kezelési útmutató GenTek™ Szögbenálló Csavarcsatornás TiBázis	49
<b>IT (Italiano)</b> Istruzioni d'uso TiBase con canale della vite angolato GenTek™	53
<b>LT (Lietuvių k)</b> Naudojimo instrukcija GenTek™ kampinio varžto kanalo TiBase pagrindas	57
<b>LV (Latviski)</b> Lietošanas rokasgrāmata GenTek™ slīpa skrūves kanāla TiBase	61
<b>NL (Nederlands)</b> Gebruikershandleiding GenTek™ TiBase met gehoekt schroefkanaal	65
<b>NO (Norsk)</b> Instruksjonshåndbok GenTek™ TiBase med vinklet skruekanal	70
<b>PL (Polski)</b> Instrukcja stosowania System łączników do kątowych kanałów na śruby GenTek™	74
<b>PT (Português)</b> Manual de Instruções GenTek™ Base de titânio de canal de parafuso angulado	78
<b>RO (Română)</b> Manual de instrucțiuni Bază din titan cu canal de șurub angular GenTek™	82
<b>SK (Slovak/Slovakian)</b> Návod na použitie Titánové bázy so skrutkovým kanálikom s nastaviteľným sklonom GenTek™	86
<b>SL (Slovenščina)</b> Navodila za uporabo Element GenTek™ TiBase za nagnjene vijalne kanale	90
<b>SV (Svenska)</b> Bruksanvisning GenTek™ vinklad skruvkanal TiBase	94
<b>TR (Türkçe)</b> Kullanım Kılavuzu GenTek™ Açılı Vida Kanallı TiBase	98

# Instruction Manual GenTek™ Angulated Screw Channel TiBase\*

\*Including final prosthetic screw

Version: 3 / 06. 2023  
Compiled at: 23. 06. 2023  
Effective date: 26. 06. 2023  
Replaces version: 2 / 11. 2022

## Indications

### Intended Use / Intended Purpose


The GenTek™ Angulated Screw Channel TiBase, an implant-level TiBase system, along with the final prosthetic screw is intended to support screw and cement-retained single-unit (engaging), multi-unit (non-engaging), and full-arch prosthetic restorations. It is used as an interface between a dental implant and an abutment restoration. The GenTek™ Angulated Screw Channel TiBase systems are intended to be used only with the corresponding GenTek™ Hexalobular Screws (final prosthetic screw and Try-in Screw).

The gluing screw is intended for laboratory use during the cementing process of the TiBase when the screw channel is straight. The gluing screw is used to block the screw channel from the potential inflow of glue while the prosthetic restoration is being glued to the TiBase, prior to final placement of the TiBase and screw into the mouth.

### Device description

Each GenTek™ Angulated Screw Channel TiBase is packaged with the final prosthetic screw and a black gluing screw. The final prosthetic screw for each system is also available as an individually packaged component. The TiBase, prosthetic screw, and gluing screw are for single use only. The TiBase and prosthetic screw are provided non-sterile and intended to be cleaned and sterilized prior to placement into the patient's mouth.

The GenTek™ Angulated Screw Channel TiBases allow the placement of the screw channel up to 25° at 360° rotation. The TiBase contains an overall post height of 7.0 mm with horizontal cutting grooves to indicate height adjustments at 4.5 mm and 5.7 mm. Additionally, the GenTek™ Angulated Screw Channel TiBases contain 3 vertical anti-rotation grooves designed to prevent the restoration from rotating. The GenTek™ Angulated Screw Channel TiBases are available with genuine connections for the following implant-level Zimmer Biomet Dental implant systems: Certain®, TSV® / Trabecular Metal™, and Ezzetic®.

For specific dental implant compatibility, please refer to the information on the individual product label and the GenTek™ components catalog.

### Materials

GenTek™ TiBases are manufactured from biocompatible titanium alloy grade 5 ELI (TiAl6V4 ELI).

## Procedure

### × Designing the Custom Restoration

The GenTek™ TiBase can be processed via a digital or traditional workflow utilizing GenTek™ Scan Bodies. The custom restoration must be constructed using dental design software. The geometry of the TiBase abutment is contained in the library file. Minimum wall thickness requirements should be controlled according to the material supplier's recommendations.

### × Milling a Custom Restoration

The custom restoration must be manufactured with milling equipment adequate for the designed geometry and material of the restoration.

### × Post-Processing the Custom Restoration

Follow the post-processing steps in accordance with the instructions issued by the material supplier.

### × Preparing the Custom Restoration for Cementing

Once the custom restoration is ready for cementing, refer to the material supplier's instructions for use for proper cementing process.

### × Preparing the GenTek™ TiBase for Cementing

Before cementing the restoration (crown or framework) to the abutment, the gluing surface of the TiBase should be thoroughly cleaned as outlined in the **Cleaning and Sterilization** section of this IFU. Blast GenTek™ TiBase in accordance with the instructions issued by the cement supplier.

The post of the GenTek™ Angulated Screw Channel TiBases has a maximum height of 7.0 mm and contains 2 cutting grooves which indicate 4.5 mm and 5.7 mm post heights. The abutment can be left at its maximum post height of 7 mm or it can be cut at one of the three indicated heights that best fits the case so that the prosthetic can be digitally designed from the corresponding library file. Abutments may be cut per standard laboratory practices and procedures and should not be cut back below the 4.5 mm groove.

### × Cementing

Follow cementing procedures in accordance with the instructions issued by the cement supplier. In order to cement the restoration to TiBase, insert the base into an implant analog and fix it in place with the gluing screw (finger-tighten). The gluing screw should be kept in place throughout the cementing process. The gluing screw should not be used when the channel is angulated on the restoration. Apply the cement in a circular fashion to the top edge of the contact surface of the TiBase and to the internal geometry of the restoration. Rotate the restoration gently as you push it over the TiBase until resistance is felt. The GenTek™ Angulated Screw Channel TiBases, single-unit and multi-unit, contain anti-rotation grooves that the restoration needs to be aligned prior to reaching its final position.

**Note: Avoid raising the restoration after the final position has been reached**, as this can compromise the integrity of the bond. The custom restoration must fit the GenTek™ TiBase without gaps and voids. Any traces of cement should be removed immediately. **The restoration must be held in position during cement removal.**

### × Polishing

Once the cement has set, remove the excess cement with rubber polishing tools. When polishing, fix the TiBase/restoration in an implant analog to prevent damaging the connection geometry.

# Instruction Manual

## GenTek™ Angulated Screw Channel TiBase\*

\*Including final prosthetic screw

Version: 3 / 06. 2023  
 Compiled at: 23. 06. 2023  
 Effective date: 26. 06. 2023  
 Replaces version: 2 / 11. 2022

### × Placing / Removing the TiBase Restoration

#### GenTek™ Angulated Screw Channel TiBase TSV® / Trabecular Metal™ and Eztetic®


**Note:** GenTek™ Hexalobular Titanium Screws (TSV® / Trabecular Metal™ and Eztetic®) are utilized for both Laboratory (Try-in) and Clinical (Final) settings for both single-unit and multi-unit restorations.


**Note:** GenTek™ Hexalobular Titanium Screws (TSV® / Trabecular Metal™ and Eztetic®) have a particular screw head design with three small external protrusions.

#### Placing in a Laboratory / Clinical Setting

1. Fully insert the Hexalobular Titanium Screw into the channel of TiBase Restoration (prior to placing into the analog/implant) using the Hexalobular Screwdriver. To ensure that the screw is properly inserted into the TiBase, screw threads must be visible at the bottom of the TiBase.  
**OPTIONAL:** If it is desired to cut the abutment down to one of the abutment heights (4.5 mm or 5.7 mm), a cutting tool can be used to cut the abutment on the machined grooves.

2. Once the Hexalobular Screw is fully inserted into the TiBase, keep pressure on the screw with the Hexalobular Screwdriver while placing the TiBase/Screw assembly onto the implant in the desired orientation. Tighten the TiBase/Screw assembly to the implant using the driver. Refer to the below table for placement torque recommendations for Hexalobular Screws.


#### GenTek™ Angulated Screw Channel TiBase Certain®


**Note:** The corresponding GenTek™ Hexalobular Titanium Try-In Screw is intended for laboratory use only. There are two different GenTek™ Hexalobular Titanium Try-In Screws (Certain®) for laboratory (Try-in) settings: one screw for single-unit restorations and one screw for multi-unit restorations. GenTek™ Hexalobular Gold-Tite® Screws are for final prosthetic use and are not intended for laboratory use.

#### Placing in a Laboratory Setting - Certain®

1. Select the proper Certain® TiBase by matching it to the implant platform diameter.
2. Place the proper Certain® TiBase into the analog in the model, line up the hex and press firmly until hearing and feeling an audible and tactile click.  
**OPTIONAL:** If it is desired to cut the abutment down to one of the abutment heights (4.5 mm or 5.7 mm), a cutting tool can be used to cut the abutment on the machined grooves.

3. Fully insert the hexalobular titanium Try-in Screw into the channel of the TiBase, using the GenTek™ Hexlobe Screwdriver and hand-tighten.

# Instruction Manual GenTek™ Angulated Screw Channel TiBase\*

\*Including final prosthetic screw

Version: 3 / 06. 2023  
Compiled at: 23. 06. 2023  
Effective date: 26. 06. 2023  
Replaces version: 2 / 11. 2022


## Placing a Clinical Setting- Certain®


**Note: The corresponding GenTek™ Hexalobular Gold-Tite® Screw is intended for final prosthetic use only. There are two different GenTek™ Hexalobular Gold-Tite® Screws (Certain®) for prosthetic settings: one screw for single-unit restorations and one screw for multi-unit restorations.**

1. Remove any provisional abutment, cover screw, healing abutment or similar from the implant upon which the TiBase is intended to be placed.
2. Confirm that the selected TiBase is the same size as the implant platform diameter.
3. Ensure that the seating surface of the implant connection is clean, free of residue, bone, or soft tissue and is dry when mounting.
4. Place the Certain® TiBase restoration onto the implant in the desired orientation, line up the hex and press firmly until hearing and feeling an audible and tactile click.

5. Once the TiBase restoration is seated into the implant, thread the Hexalobular Gold-Tite Screw into the TiBase restoration, tighten the TiBase restoration/screw assembly to its prosthetic torque using the hexalobe screwdriver. Refer to the below table for placement torque recommendations for Hexalobular Screws.

## × Removing the TiBase Restoration

1. Should the TiBase/Screw (TSV®/Trabecular Metal™ and Eztetic®) system require disassembly from the implant, an abutment removal tool is not required. The additional feature on the Hexalobular Titanium Screw is intended to assist in the removal of the TiBase during standard screw removal.

**Note: Should the restoration remain assembled to the implant once the screw has been loosened, direct manual manipulation of the restoration with or without a tool such as a hemostat, may be used to facilitate full removal.**

Placement torque recommendations for GenTek™ Hexalobular (laboratory Try-in or final prosthetic) Screws:

**Note: Only corresponding GenTek™ Hexalobular Screws for each connection (Certain®, TSV®/Trabecular Metal™, and Eztetic®) must be used with the corresponding GenTek™ Angulated Screw Channel TiBase. Hex drive screws from other GenTek™ TiBase systems must not be used with the GenTek™ Angulated Screw Channel TiBase system.**

Connection*	Angulation Range	Laboratory Try-in Screw	Final Prosthetic Screw Torque	GenTek™ Hexalobular DriverScrew*
Certain®	0 – 25 degrees	Hand-tighten**	20 Ncm	ZFX02HLD21, ZFX02HLD28, ZFX02HLD36
TSV® / Trabecular Metal™	0 – 25 degrees	Hand-tighten**	30 Ncm	ZFX02HLD21, ZFX02HLD28, ZFX02HLD36
Eztetic®	0 – 25 degrees	Hand-tighten**	30 Ncm	ZFX02HLD21, ZFX02HLD28, ZFX02HLD36

\*Applicable to single-unit and multi-unit for each connection

\*\*Minimum of 5 Ncm; Maximum of 15 Ncm

## Cleaning and Sterilization

The GenTek™ Angulated Screw Channel TiBase is non-sterile. Prior to placement of the GenTek™ TiBase/restoration and prosthetic screw into the implant, the GenTek™ TiBase assembly must first be cleaned and sterilized.

### Step 1: Manual Cleaning

1. Rinse the instruments under running cold tap water for 2 minutes.
2. Hand wash instruments in a bath with demineralized water for 5 minutes. Remove debris with a soft brush. Cavities should be rinsed with a syringe. This step should be repeated until all debris is removed.
3. Clean the instrument in an ultrasonic bath using a mildly alkaline cleaner (1.5%) at a frequency of at least 40 kHz for 20 minutes.
4. Rinse thoroughly with demineralized water for 1 min.

### Step 2: Sterilization

One of the following sterilization methods is recommended to be used:

#### × Method 1:

Full cycle pre-vacuum steam sterilization at a temperature of 134 °C (273.2 °F) for 3 minutes, 4 pulse minimum.

# Instruction Manual GenTek™ Angulated Screw Channel TiBase\*

\*Including final prosthetic screw

Version: 3 / 06. 2023  
Compiled at: 23. 06. 2023  
Effective date: 26. 06. 2023  
Replaces version: 2 / 11. 2022

- × Method 2: (alternative to method 1)  
Full cycle gravity steam sterilization at 134°C (273.2°F) for 3 minutes minimum.

### Step 3: Drying

After sterilization treatment it is recommended to dry the products at a temperature of maximum 134°C (273.2°F) for 10 minutes. The products are ready for use 5 minutes after completion of the drying process.

**Note: Follow the loading instructions issued by the sterilizer manufacturer.**

### Storage and Handling

The GenTek™ TiBase should be stored at room temperature.

### Contraindications

The custom abutment and/or restoration that is designed shall not be angulated more than 20° off the implant axis.

### Known Side Effects

Allergies to the alloy or contents of the alloy may very rarely occur.

### Warnings

- × The GenTek™ TiBase must only be used and handled by dental professionals who have approved training in dentistry.
- × GenTek™ TiBase may only be used for Zimmer Biomet Dental implant systems.
- × The black gluing screw is not intended to be used in the patient's mouth.
- × The black gluing screw should not be used when the screw access on the crown is at an angle.
- × The GenTek™ TiBase and prosthetic screw are small and, therefore, they must be handled with caution to avoid swallowing or inhaling by the patient.
- × The use of an incorrect prosthetic screw may result in damage to the GenTek™ TiBase, the prosthetic screw and / or the implant.
- × The use of a prosthetic screw torque higher than the recommended value on the label may result in damage to the GenTek™ TiBase driver, the prosthetic screw and / or the implant. The use of torque values lower than those recommended may result in loosening of the GenTek™ TiBase and / or the prosthetic screw.
- × The screwdriver must be in good condition . A worn screwdriver can lead to over-rotation of the hex in the screw head which can result in the screw not being fully tightened and / or not being able to be removed.
- × For any lab work and any try-in of the GenTek™ TiBase, a lab screw (non-Gold-Tite) must be used. The included prosthetic screw may only be used to permanently fix the finished restoration.

### Explanation of symbols

Non-sterile

Do not use if package is damaged

Manufacturers Information Address of legal manufacturers facility

Don't use twice

Recommended torque to be applied for final placement of the included screw

Recommended to fix the Try-In Screw hand-tight

LOT-Number, Symbol followed by Production Date of the lot

Medical Device

Reference Number, Symbol followed by Manufacturer's Item Number

Quantity / Unit of package content

Unique Device Identifier

Consult instructions for use [www.zfx-dental.com](http://www.zfx-dental.com)

Consult instructions for use for mating hexalobular screws and angulated screw channel TiBase abutments.

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Zfx GmbH