


ProGraft system

CE
0123

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ProGraft system allows effective infusion and hydration of particulate bone substitutes and easy direct grafting in bone defects in oral surgery.

The system comprises a device with a blue piston to graft the bone substitute, a funnel, a device with a green piston to draw the sterile saline solution, a connector to connect the devices, and a yellow piston to prepare the grafting device for reuse.

- The ProGraft system device has been tested with BIO-OSS hydroxyapatite. **Before use of different particulate bone substitutes, it is essential to conduct a functional test.**
- Before use, check that the package is intact.
- Do not use the product if the package is damaged.
- The ProGraft system device must be used exclusively by skilled medical staff.
- Always wear sterile gloves while handling the ProGraft system device and strictly follow the procedures in order to guarantee sterility of the biomaterial treated.
- The ProGraft system device is single-use and absolutely must not be re-sterilized.
- **Warning:** if it is reused on another patient, there is a risk of cross-contamination and loss of the performance and functional characteristics of the device.
- Discard after use in special medical waste containers in compliance with the regulations in force.
- **META does not shoulder any responsibility for improper use of the product.**

- 1 Open the package on the side indicated and deposit the ProGraft system device on the sterile field always handling it in absolutely sterile conditions.
- 2 Grip the device with the green piston and draw the sterile saline solution pulling the piston backward but not fully out (Figure 1).

- 3 Grip and position the funnel on the end of the device with the blue piston, check that the piston is in backward position with the shiny notch visible, and pour in the particulate bone substitute helping it along by lightly tapping on the funnel making sure that the device is completely filled (Figures 2A, 2B and 2C).

- 4 Fit the connector in place of the funnel and join the device with the blue piston to the device with the green piston containing the sterile saline solution (Figures 3A and 3B).

- 5 For better hydration of the bone substitute, position the two devices vertically with the blue ProGraft device at the top. Gently push the green piston until reaching the base of the device so as to hydrate the bone substitute in the bone graft device with the blue piston (Figure 4).

- 6 Complete hydration of the bone substitute with sterile saline solution is visually indicated thanks to the transparent body of the device with the blue piston. Remove the connector from the device with the blue piston, let the excess liquid flow out and continue with extrusion of the bone substitute directly in the surgical site (Figures 5A and 5B).

- 7 For repeated use of the bone graft device, extract the blue piston, insert the yellow piston provided and perform a complete cycle to clean the device. Reposition the blue piston and follow the procedure in point 2.

Optimal functionality of the bone graft device with the blue piston is guaranteed for 3 uses within the same surgical operation.

Warnings:

The surgeon must establish whether the patient is fit for bone regeneration operations.

Only breaching the contents of this document can you avoid errors and ensure proper functioning of our product.

These instructions for use together with further information on ProGraft system are available on the website www.metahosp.com

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