

Identium[®] Heavy



Here used in combination with Identium[®] Light in the double mix impression technique

Identium® – the best of two worlds.

Identium[®] Heavy is a high viscosity, addition-curing elastomeric impression material based on Vinylsiloxanether[®]

Identium[®] Heavy has been specially developed for the one-step impression technique and combines the best of two worlds: all the advantages of polyether materials coupled with those of A-silicones.

The highly elastic properties of Identium[®] Heavy ensure dimensionally accurate recovery and very easy removal. Identium[®] Heavy is also odorless and has a neutral taste.

Based on the entirely new material of Vinylsiloxanether[®], Identium[®] Heavy provides the following benefits:

 Excellent flowability combined with exceptional hydrophilicity (lowest achievable contact angle of less than 10° after 1 second), ensuring optimal flow behavior in a moist environment.

Intro pack31 x 300 ml base, 1 x 62 ml cata1 x 50 ml Light body, 6 mixing	REF 28237					
Refill pack72 x 300 ml base, 2 x 62 ml cata	24 ml _{alyst}	REF 13776				
Economy pack 2 6 x 300 ml base, 6 x 62 ml cat		REF 13777				
Current prices are available from our distributors. Not all products/packages are available in all markets.						

Vinylsiloxanether[®] (VSXE) Precision impression material

SUGGESTED IMPRESSION TECHNIQUE

One-step impression technique (simultaneous technique)

INDICATIONS

Functional impressions

Implant impressions

Impressions for crowns/bridges/inlays/ onlays and veneer preparations

Impressions about full or partial denture

highly recommended recommended



Precision impression material (VinyIsiloxanether®) The Identium® Heavy Advantages and Benefits

Product Advantages

- Excellent flowability
- Exceptional hydrophilicity
- Optimal elastomeric properties
- Balanced setting behavior
- Odorless and neutral taste
- Accurate model fabrication


Contact angle measurement

Easy handling

Benefits to the Practitioner

- Material flows into the narrowest crevices with high positional stability, constantly guaranteed throughout the entire working time.
- Optimal flow behavior in a moist environment, lowest achievable contact angle (less than 10° after 1 second!).
- Dimensionally accurate recovery and easy removal from the patient's mouth.
- Long working time with short intraoral setting time (double-snap effect).
- More comfortable for the patient: no gag reflex, no erratic movement.
- Minimum risk of fracture when removing the impression.
- In the practical Plug & Press[®] system; compatible with commercially available disinfectants (no swelling), can be poured immediately.

Identium[®]. A perfectly balanced impression material.



Identium® after 1 sec.: contact angle < 10°

Identium[®] Heavy is a high viscosity impression material with high final hardness available in foil bags. It is particularly suitable for the double mix technique (in combination with low viscosity Identium[®] Light), for taking impressions for the fabrication of crowns/bridges, inlays/onlays and veneers, for functional and implant impressions as well as impressions for removable restorations.

Technical Data

Product name	Mixing ratio	time at 23° C/	Working time at 35° C/ 95° F ≤	Intraoral setting time at 35° C/ 95° F ≥	Total setting time* ≥	Hardness (approx.) Shore	dimensio- nal change (max.) %	Elastic Recovery test (approx.) %	Strain in compres- sion (approx.) %
Identium® Heavy	5:1	2 minutes	80 seconds	3 minutes 30 seconds	5 minutes 30 seconds	A 60	- 0.20	99.0	2.8
Product name	ISO 4823	Consistency (approx.) mm	Mixing ration content vol unit		Mixing techniq	ue	Mixing ele- ment	Color	Odor/taste
Identium® Heavy	Type 1, heavy-bodied	32	5:1, 362 ml foil bag		Plug & Press [®] dispenser or other automatic mixing and dispensing unit		Dynamic mixer	Light blue	Neutral

* Total setting time (removal from the mouth) from the start of mixing.