INSTRUCTIONS FOR USE

DESCRIPTION

D.T. LIGHT-POST® **ILLUSION™ X-RO®** is a yellow translucent fiber post with a colored resin matrix which becomes translucent in the mouth. The post color will re-appear on command under cold air, water spray if re-treatment is required.

The radiopacity level of the **D.T. LIGHT-POST® ILLUSION™ X-RO®** is 340% Al eqv. (ISO 4049)

	■ 0,5	■ 1	■ 2	■ 3
∅ apical tip	0.80	0.90	1.00	1.20
ØPost head	1.25	1.50	1.80	2.20
Double taper	.02/.04	.02/.06	.02/.08	.02/.10

INDICATIONS FOR USE

In case of insufficient residual tooth structure (<4 mm), the post is needed to support and secure the coronal restoration.

CONTRAINDICATIONS

If coronal tooth structure is less than 1.5 mm.

PRECAUTIONS

The post must be cleaned with alcohol before insertion in the canal. Avoid touching the posts with your fingers after cleaning.

Shortening the post should be done outside of the mouth.

The use of rubber dam is recommended.

The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.

ADVERSE REACTIONS

None known to date.

STEP-BY-STEP INSTRUCTIONS

- Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
- Select the correct **D.T. LIGHT-POST® ILLUSION™ X-RO®** size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
- Determine the length of the coronal build-up, which will account for 1/5 of the final post length. Remove the root filling with the **D.T. Universal Drill** (rotation speed 1000-2000 rpm) as deep as necessary for the post to be inserted to 2/5 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
- Shape the canal with the **D.T. finishing Drill** corresponding to the selected **D.T. LIGHT-POST® ILLUSION™ X-RO®** post size, (rotation speed 1000-2000 rpm).
- Try in the post, to check that the post is properly seated in the canal
- Shorten the post to its final length with a diamond disc outside of the mouth. **Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.**
- Clean the post with alcohol.
- Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. **Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.**
- Apply 2 coats* of the primer with a thin brush (e.g. Composibrush RTD) in the post space. Remove excess using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
- Apply a single coat of primer* to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth.
- 11a.** If using a dual cure resin cement for both cementation and core build up such as Corecem™ (recommended technique): Mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post immediately. Light-cure for 40 to 60 seconds. Continue applying the core build up using a plastic Coreform™ matrix, or directly on the tooth.
- 11b.** If using separate cement and a core build up resins : Mix

the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).

- 11c.i.** If using a self etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.
- 11c.ii.** After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).
- Directly model the build-up using a composite core build-up material*. R.T.D. recommends that the coronal end of the post be covered by core composite. In cases where a prosthetic crown is not indicated, a minimum of 1 mm of composite restorative should cover the post and that should be an appropriately-shaded, polishable, aesthetic composite. If the **ILLUSION** post cannot be adequately covered with composite, then a non **ILLUSION** version can be utilized.

OTHER RECOMMENDATIONS

A thick rubber dam can momentarily prevent the post to change color.

D.T. LIGHT-POST® ILLUSION™ X-RO® are single use devices. They must not be re-used to avoid the risk of contamination.

D.T. LIGHT-POST® ILLUSION™ X-RO® can be cleaned with alcohol. They can neither be sterilised nor disinfected in thermodesinfector.

The use of autoclave is not recommended.

A re-access carbide kit is available if needed. Contact your RTD distributor.

D.T. Universal Drill and D.T. Finishing Drill: these instruments are supplied non-sterile. Disinfect and sterilize instruments before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermodesinfector. Inspect instruments visually before each use: change to new drills when they are damaged, corroded or dull. The Starter drill may be used only to remove root canal filling partially, and the finishing drills may be used only to prepare root canal.

Starter Drill and Finishing drills must be sterilized individually in an autoclave with the following settings**:

- Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use (e.g EN ISO 11607-1).
- Autoclave: type B complying with regulation standard in force (e.g EN 13060).
- Sterilization temperature: 134°C - sterilization time: 18 minutes.

Use a new drill after 12 -15 uses.
For profoundly flared or ovoid canals, the use of accessory post, such as Fibercone™, is indicated for best results.

STORAGE

Store the sterilized components in a dry, dust-free place. If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch or re-sterilize according to the protocol defined before.

Do not store the post in direct sun light or artificial light.

SAFETY AND LIABILITY

RTD declines any liability or compensation for possible damage due to:

- use of instruments which do not belong to the system and which may impair the function,
- failure to observe the instructions for use.

The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the instructions for use. Selection of the correct drill and post is the responsibility of the user.

For dental use only.

* Refer to the manufacturer's instructions for use.
** Please refer to your autoclave manufacturer for specific sterilization instructions if needed

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